


Nanodyne Replacement Illuminator for Zeiss Axiovert Microscope: Included Items



**Additional Items Included But Not Shown:**

PN 10490 Hex Key 5/64 inch for set screw



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	4/3/2020
PN 10760 Zeiss Axiovert Illuminator System	REV 3
SHEET 1	OF 1