

**Nanodyne Replacement Illuminator for Zeiss Axiostar Microscope: Included Items**

**Additional Items Included But Not Shown:**

PN 10457 Hex Key 2.5mm Short Arm

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



Tape to secure wires.

Han-Seung Yang 4/4/2020	
PN 11718 Zeiss Axiostar illuminator System	REV 1
SHEET 1 OF 1	

