

Nanodyne Replacement Illuminator for Zeiss AXIO Vert.A1 Microscope Illuminator: Included Items

Additional Items Included But Not Shown:

- PN 10490 Hex Key 5/64 inch for set screw
- PN 10747 Hex Key 3mm



PN 10736 Rubber plug to block unused AC power receptacle.



PN 10733 Power Supply and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



Tape to secure wires.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang		6/19/2020	
PN 11773 Zeiss AXIO Vert.A1 illuminator System		REV 1	
SHEET 1		OF	1