


Nanodyne Replacement Illuminator for Zeiss Axiovert 25 Microscope: Included Items



PN 10983 T2PI Main
Assy R3 with Lens

PN 10467 Brass Tip
Set Screw 8-32

PN 10756 Zeiss Axiovert Mounting Ring

Additional Items Included But Not Shown:

PN 10490 Hex Key 5/64 inch for set screw



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	4/2/2020	
PN 11935 Zeiss Axiovert 25 Illuminator System		REV 3
SHEET 1 OF 1		