


Nanodyne Replacement Illuminator for Zeiss Axioskop 2 Microscope: Included Items



Additional Items Included But Not Shown:

PN 10490 Hex Key 5/64 inch for set screw



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang	6/25/2021
PN 12066 Zeiss Axioskop 2 Illuminator System (low power)	REV 1
SHEET 1 OF 1	