

Nanodyne Replacement Illuminator for Zeiss Axioskop 2 Microscope: Included Items

PN 11145 T2PIHP3 Illuminator Assy w Offset Lens

PN 11243 Knob 0.236 in shaft
0.625 dia CLR Gloss


PN 11170 Pot-Cable Assy
(detabbed, 5mm shaft cut)
(15 in cable)

PN 10467 Brass Tip
Set Screw 8-32 x 0.188

PN 11247 Zeiss HAL 100 Adapter

Additional Items Included But Not Shown:

PN 10490 Hex Key 5/64 inch for set screw



PN 11163 Power Supply 5V 2.1A and PN 10734 Power Cable 1.35mm ID/3.5mm OD x USB A 6 ft.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang	6/25/2021
PN 12067 Zeiss Axioskop 2 Illuminator System (high power)	REV 1
SHEET 1 OF 1	