

Nanodyne Replacement Illuminator for Olympus BX40 Microscope - Included Items



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output



Tape to secure wires.



PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.




PN 10736 Rubber plug to block unused AC power receptacle.

Power supplies subject to substitution without notice due to availability issues and changes in regulations.



PN 10908 BX40 Pot Plate - Cable Assy (9 inch long cable)

PN 10334 Olympus BX40 Series Illuminator Assy



Wayne Bonin	1/9/2019	
PN 10450 Olympus BX40 Full Illuminator System		REV 3
SHEET 1 OF 1		