


Nanodyne Replacement Illuminator for Olympus BX41/45 Microscope - Included Items

PN 10334 Olympus T2PI
BX40 Series Illuminator Assy



***Note: Nanodyne's Illuminator System PN 10472 is used for both the Olympus BX41 and BX45.**

PN 10603 BX41 Pot-Cable Assy

Pot Knob from BX41
(save from original microscope)

Additional Items Included But Not Shown:

PN 10456 Hex Key 1.5mm
3 pcs PN 10610 cable ties
3M VHB tape

Not Included:

Knob - Save from original BX41 microscope

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot



Tape to secure wires.



PN 10736 Rubber plug to block unused AC power recepticle.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	2/3/2019	
PN 10472 Olympus BX41-45 Illuminator system		REV 6
	SHEET 1	OF 1