


**Nanodyne Replacement Illuminator for Olympus BH-2 Microscope - Included Items**



**Additional Items Included But Not Shown:**  
 PN 10457 Hex Key 2.5mm (for M3 mounting screws)


**Not Included - Save From Original:**  
 Knob from old slide pot



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output



PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A 6 foot



PN 10736 Rubber plug to block unused AC power receptacle.

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin	1/8/2019	
PN 10474 Olympus BH-2 Full Illuminator System		REV 5
	SHEET 1 OF 1	