


Nanodyne Replacement Illuminator for Nikon 50i Microscope - Included Items



Additional Items Included But Not Shown:

- PN 10456 Hex Key 1.5mm (for pot knob)
- PN 10566 Hex Key 2.5mm long arm ball end (for M3 mounting screws in original illuminator collector lens assembly)
- PN 10490 Hex Key 5/64 inch (for 2-56 mounting screws)

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.



PN 10733 Power Supply - XP Power 5V 1A and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.



Tape to secure wires.



PN 10736 Rubber plug to block unused AC power recepticle.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.

Wayne Bonin		2/13/2019
PN 10558 Nikon 50i Illuminator System		REV 6
SHEET 1		OF 1