



Nanodyne Replacement Illuminator for Olympus BH Microscope - Included Items

Additional Items Included But Not Shown:
 PN 10457 Hex Key 2.5mm (for M3 mounting screws)

Not Included:
 Intensity Adjust Knob - Save From Old Slide Pot



PN 11002 Pot Cable Assy

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies are subject to substitution without notice due to availability issues and changes in regulations.



PN 10733 Power Supply - XP Power 5V 1A Fixed Blade US Input/USB Output and PN 10734 Cable Assy 1.35mm ID x 3.5mm OD RA plug to USB A, 6 foot.

Wayne Bonin	1/5/2018
PN 10716 Olympus BH Full Illuminator System	REV 5
SHEET 1	OF 1