



PN 11145 T2PIHP3 Illuminator Assy
w Offset Lens

PN 10756 Zeiss Axiovert Mounting Ring
modified for Leitz Wetzlar SM-LUX-POL

OEM Mounting plate
(save from original microscope)



PN 11028 Power Supply - CUI 5V 2.5A 2.1 x 5.5mm RA plug



PN 10736 Rubber plug to block unused AC power recepticle.



Tape to secure wires.

***Note that the illuminator mounted adjustment potentiometer is standard on this unit. Ask about a remote adjustment pot if desired.**

Additional Items Included But Not Shown:
 PN 10709 Hex Key 3/32 Short Arm (for 4-40 mounting screws)
 PN 10337 Hex Key 2mm Short Arm (for set screw)

The illuminator may be powered by plugging the cable into the power supply provided, or into a suitable USB port on a computer or other device.

Power supplies subject to substitution without notice due to availability issues and changes in regulations.

Han-Seung Yang		8/27/2018
PN 11220 Leitz Wetzlar SM-LUX-POL 100W illuminator System		REV 1
SHEET 1	OF	1